



The Pellevé™ Wrinkle Reduction System

- + Immediate, measurable improvement
- + No downtime
- + Virtually no pain
- + Safe and simple to use
- + Innovative
- + Flexible and versatile

It All Adds Up

pellevé™
Elegance through science.

ellman®
Experts in Precision Surgery.

The Pellevé™ Wrinkle

Noticeable Results

The Pellevé™ Wrinkle Reduction System is an innovative, non-invasive skin tightening treatment for mild to moderate facial wrinkles. Pellevé™ uses advanced radiowave technology to induce collagen contraction and synthesis without damaging the epidermis, so your patients can look and feel refreshed – naturally.


Immediate, measurable improvement

- Visible and lasting effects in nearly all patients
- Better skin quality and reduction in skin laxity

Safe and effective

- Non-invasive
- No anesthetic required
- No skin cooling needed
- Minimal discomfort and side effects

Sustained results over time¹



* Response rate of total study population after one treatment with Pellevé™ (n=93)¹



Before | After (Immediately)

Notice elevation of ptotic brow, reduced crow's feet and smoother under-eye area.

Photo courtesy of Dr. D. Richey



Before | After (180 Days)

Notice improved texture of the cheek, reduced bag under the eye, and reduced lines around the mouth.

Photo courtesy of Dr. A. Rusciani

1. Rusciani A, Curinga G, Menichini G, Alfano C, Rusciani L. Nonsurgical tightening of skin laxity: A new radiofrequency approach. *J Drugs Dermatol*. 2007;6(4):377-382.

e Reduction System

+ Flexible System

NEW Pellevé™ S5 Generator – complete aesthetic versatility with the proven Surgitron® platform

The new Pellevé™ S5 provides maximum control in precision cutting and energy delivery, plus more versatility than other energy-based technologies², including:

- **Simple** – plug in the GlideSafe™ handpiece, press the button and begin
- **Quick** – handpiece is ready instantly, no waiting for unit to heat up
- **Versatile** – a comprehensive aesthetics platform with hundreds of electrodes/accessories available for a variety of surgical procedures
- **Flexible** – ideal for use in both office and clinic environments



“The results Pellevé offers for the eyes I cannot get with any other device – nothing compares to it for the eyelid area because the tip is finer and more delicate.”

– Dr. Flor Mayoral, dermatologist, Miami, USA

Full surgical system with S5 Generator, Surg-e-Vac™ smoke evacuation unit and Cart.




The Pellevé™ Wrinkle

Innovative Technology – the new Glide

It's all under your control – the **NEW** GlideSafe™ Fixed-Life Reusable Handpiece

GlideSafe™, the newest innovation from the Pellevé™ brand, provides continuous energy flow, not pulsed. It features an **innovative sensor** that will cut the energy flow when the device senses loss of contact with the skin – virtually eliminating the risk of shocking – so it's **more comfortable than ever** for your patients.



“Fixed-Life Reusable” defined

- Provides **3 cumulative hours of activation time** – treatments can be non-contiguous
- Handpiece can be used for multiple procedures and/or patients – great for touch-ups
- Running time can be started and stopped as desired – treat **how** and **when** you want
- The 3-hour GlideSafe™ LifeTime is consumed only while RF energy is flowing



eReduction System

eSafe™ Handpiece*



Proprietary
gold-plated
tip design



Elegant design, flexible use

- All faces are different – now you can vary each treatment as needed
- Compact and maneuverable – easy to hold
- Simple to use – fingerswitch or footswitch operation
- Great for sculpting – easily gets into the contours of the face

Available in 4 sizes for added flexibility

- 7.5mm for tightly confined areas
- 10mm for less confined areas
- 15mm for large open areas
- 20mm for larger open areas



“The tip is perfect for the eyelid, under-eye or the fine lines on the top lip. Results are very quick with Pellevé and patients love that!”







– Dr. Rita Rakus, aesthetic doctor, London, UK

* Patent pending

The Pellevé™ Wrinkle

Four Sizes, Multiple Uses

GlideSafe™ Handpiece Uses

		7.5mm	10mm	15mm	20mm	
Surface Area	Less					More
Skin Contact	Less					More
Procedure Time	Slower					Faster
Patient Turn-Over	Less					More
Uniformity of Treatment	Less					More
Applications	Treat Tightly Confined Areas					Treat Larger Unconfined Areas
Ideal Usage Areas		Periorbital Glabellar	Periorbital Glabellar Perioral	Forehead Cheeks Jawline Perioral	Forehead Cheeks Jawline	

Pellevé™ Accessories

Fluke 62 Mini Infrared Thermometer

TEMP62

Pellevé™ Treatment Gel

2oz DHPG2

8oz DHPG8

Pellevé™ Cold Pack

PPK01



e Reduction System

+ Great Benefits

Profitability for your practice


- Add a new cosmetic technology to your business quickly and affordably
- Generate new revenue with the wrinkle-reducing treatment patients want
- Achieve rapid return-on-investment with one system that supports a range of different surgical procedures

Become a Pellevé™ Partner

- Achieve optimal results with an exclusive education and training program
- Utilize the marketing tools available to you at pelleve.com
- **Find a Physician** listing at pelleve.com
- Ongoing support from Ellman International

It Certainly Adds Up...

- + **Efficacy** – 87% of patients experienced fewer overall facial wrinkles at 6 months¹
- + **Patient comfort** – virtually no downtime or pain
- + **Safety and innovation** – new GlideSafe™ handpiece features an innovative sensor to virtually eliminate the risk of shocking your patient
- + **Simplicity** – it's literally “plug and play”
- + **Flexibility and versatility** – a comprehensive aesthetics platform with hundreds of electrodes/accessories available for a variety of surgical procedures



GlideSafe™ 15mm
handpiece

To learn more about the **Pellevé™ Wrinkle Reduction System**, visit pelleve.com

pellevé™
Elegance through science.

1. Rusciani A, Curinga G, Menichini G, Alfano C, Rusciani L. Nonsurgical tightening of skin laxity: A new radiofrequency approach. *J Drugs Dermatol*. 2007;6(4):377-382.


Complete Aesthetics Platform

The new **Pellevé™ S5** is a comprehensive aesthetics platform with hundreds of electrodes/accessories available for a variety of surgical procedures.

Pictured: Full surgical system with S5 Generator, Surg-e-Vac™ smoke evacuation unit and Cart (20" wide x 20" deep x 40" high).

"In addition to treating wrinkles, I use the S5 for surgical cases in my office – it's great for telangiectasias – and the ability to excise or biopsy lesions on the skin and mucosal surfaces, without suture repair, has saved me time and money without compromising patient outcome."

– Dr. Michael Stampar, facial plastic surgeon, Punta Gorda, USA



pellevé™
Elegance through science.

Important Safety Information

Pellevé™ is indicated for the non-ablative treatment of mild to moderate facial wrinkles for skin types I–IV. Possible side effects of Pellevé™ are usually treatment-site related and include mild discomfort during the procedure localized within the treatment area. Mild swelling (edema) and redness may occur. Initial studies indicate that these typically resolve within 2 to 24 hours.

Please see Instruction Manual for more information.